



PURPOSE

Cytometry is a powerful analytical tool for the analysis of multiple biological parameters of individual immune cells within heterogeneous cell populations [1,2]. The implementation of standardized cytometry-based biomarker assays in clinical trial support remains a challenge due to inter-instrument variations and the limited stability of clinical specimens.

Chipcytometry is an image cytometry platform enabling long-term sample storage. Here we present validation data (precision, stability) on three assays (regulatory T-cells (Treg), T-helper 17 cells (Th17) and plasma cell assays), following industry guidelines using an automated Chipcytometry system (CYTOBOT™).


METHODS

- Three assay panels were developed using commercially available, directly labelled antibodies against 7 markers (see table below).
- Human PBMC samples from 5 healthy donors were loaded in microfluidic cassettes (ZellSafe™) for precision and stability assessments.
- Assays were validated for Treg, Th17 and plasma cells, according to industry guidelines [3,4].
- Data processing of images was performed using an automated gating algorithm (ZellScanner™ software) and the statistical programming language R [5].

Cell Type	Markers
Plasma Cells	CD19/CD27/CD38
Th17 Cells	CD3/CD4/CD161
Treg Cells	CD3/CD4/Foxp3



PRINCIPLE



INSTRUMENT



RESULTS - PRECISION

Determination of the intra-assay precision showed a mean CV of 14.6% (plasma cell assay), 16.2% (Th17 assay) and 19.1% (Treg assay). The mean inter-assay CV, which comprises both within- and between-batch variability was 19.3% (plasma cell assay), 29.6% (Th17 assay) and 22.5% (Treg assay).


Parameter	Assay	Donor				
		A	B	C	D	E
Intra-Assay Precision						
Plasma Cell Assay	1	15.7	18.3	12.0	9.8	14.2
	2	18.3	18.3	10.0	ND	ND
Th17 Assay	1	12.8	13.3	15.5	10.2	8.2
	2	24.0	29.6	15.9	ND	ND
Treg Assay	1	26.7	32.3	19.5	13.2	12.6
	2	39.6	3.9	4.9	ND	ND
Inter-Assay Precision						
Plasma Cell Assay		25.1	22.6	10.1	ND	ND
Th17 Assay		34.5	27.8	26.6	ND	ND
Treg Assay		30.2	20.7	16.6	ND	ND

RESULTS - STABILITY

Stability analysis of 5 PBMC samples from 5 donors, stored in ZellSafe™ cassettes at 4° C for 3 months, revealed that 5/5 samples (Treg assay) and 4/5 samples (plasma cell assay and Th17 assay) were within the 25% range from baseline level, which was within predefined acceptance criteria.

Staining on PBMCs from a single preparation (baseline vs. 3 months)

Stability T = 0



Stability T = 3 months

CONCLUSIONS

- Three Chipcytometry assays for quantification of Treg, Th17 and plasma cells in human PBMC samples were deemed fit-for-purpose after validation.
- Chipcytometry is an alternative cytometric platform, which allows banking of PBMC samples, permitting retrospective multi-parameter batch analysis for novel biomarkers that may not have been envisioned at the beginning of a clinical trial.
- Further characterization of this method will determine its value for immunomonitoring purposes in the clinical development of immune-based therapeutics.

REFERENCES

- [1] Hennig et al. Cytometry A. 2009: A versatile platform for comprehensive chip-based explorative cytometry.
- [2] Happel et al. Science Translational Medicine. 2014: Pulmonary transplantation of macrophage progenitors as effective and longlasting therapy for hereditary pulmonary alveolar proteinosis.
- [3] O'Hara et al. Journal of Immunological Methods (2010): Recommendations for the validation of flow cytometric testing during drug development: II assays.
- [4] Wood et al. Clinical Cytometry, Part B 2013: Validation of Cell-based Fluorescence Assays: Practice Guidelines from the ICSH and ICCS.
- [5] R Core Team (2016). R: A language and environment for statistical computing. R Foundation for Statistical Computing, Vienna, Austria. URL <https://www.R-project.org/>.