



PURPOSE

Mesenchymal stromal (stem) cells (MSCs) present promising opportunities for tissue engineering and regenerative medicine. They can be isolated from various tissues, e.g. bone marrow, adipose tissue, placenta etc. Here, we report for the first time in situ detection of cells expressing this marker combination in human placenta using high-parameter immunohistology.

METHODS

Placenta samples were embedded in Tissue-Tek® and cryopreserved. All chips were shipped and stored at 4 degree C. The applied antibodies were CD105-PE (clone 43A3, Biolegend), CD73-PE (AD2, Biolegend), CD90-PE (5E10, Biolegend), CD34-PE (561, Biolegend), CD19-PE (HIB19, eBiosciences), CD14-PE (RMO52, Beckman Coulter), CD45-Alexa Fluor 488 (HI30, Biolegend) and HLA-DR-BUV395 (G46-6, BD Biosciences). Antibodies were tested with purified BM-MSCs. Cell segmentation, fluorescence value calculation and data analysis of all chips were performed with Chipcytometry software.

RESULTS I





Fig. 1 Staining of the whole human placenta cryosection (5 μm) with CD90 (left) and expression of the entire marker set by three MSC areas (right).

Acknowledgement:



PRINCIPLE



INSTRUMENT



CHIPCYTOMETRY 95-PLEX ASSAY

| | | | | | | | |
|-------|-------|--------|-------|----------------|----------------|----------|-------------|
| CD3 | CD21 | CD40 | CD73 | CD152 (CTLA-4) | CD206 | Granz. B | IL17A |
| CD4 | CD24 | CD45 | CD80 | CD161 | CD244 | Helios | IL17F |
| CD5 | CD25 | CD45RA | CD81 | CD163 | CD257 (BAFF) | HLA-DR | Ki-67 |
| CD8 | CD27 | CD45RO | CD86 | CD172a/b | CD273 (PD-L2) | IFNγ | LC (κ) |
| CD10 | CD28 | CD54 | CD90 | CD183 (CXCR3)* | CD274 (PD-L1) | IgA | LC (λ) |
| CD11b | CD29 | CD56 | CD95 | CD184 (CXCR4)* | CD278 (ICOS) | IgD | pan Cytok. |
| CD11c | CD30 | CD57 | CD105 | CD185 (CXCR5)* | CD279 (PD-1) | IgG | RORγ(t) |
| CD14 | CD31 | CD64 | CD115 | CD193 (CCR3)* | CD319 (CRACC) | IgM | T-bet |
| CD15 | CD32* | CD66b | CD123 | CD194 (CCR4)* | CD326 (EpCAM) | IL1b | TNFα |
| CD16 | CD34 | CD68 | CD127 | CD195 (CCR5)* | AIOLOS (IKZF3) | IL8 | Vimentin |
| CD19 | CD38 | CD69* | CD138 | CD196 (CCR6)* | FoxP3 | IL10 | Zap-70 |
| CD20 | CD39 | CD71 | CD141 | CD197 (CCR7)* | GM-CSF | IL12 | *live stain |

RESULTS II

Human placenta cryosections deriving from early terminations (up to 12th week of pregnancy) as well as from full term pregnancies were analysed. According to the classical ISCT definition, we found six MSCs on tissue section M126694, three MSCs on M126688 and no MSCs on M126696 (Tab 1). One representative staining of section M126688 is depicted in Figure 1. Since MSCs might not express all typical surface molecules CD105, CD73 and CD90 in vivo, we also searched for single- and double-positive cells (Tab 1). Between one and seven CD73+ CD90+ regions could be detected on the different sections, while none or up to two CD73+ CD90+ areas, could be localized, respectively.

RESULTS III

| Chip ID | Week of gestation | CD73+ CD105+ CD90+ lineage - | CD73+ CD90+ lineage - | CD73+ CD105+ lineage - | CD73+ lineage - | CD105+ lineage - |
|---------|-------------------|------------------------------|-----------------------|------------------------|-----------------|------------------|
| M126688 | 11th | 3 | 3 | 2 | 3 | 5 |
| M126694 | 39+0 | 6 | 7 | 0 | 0 | 0 |
| M126696 | 39+6 | 0 | 1 | 1 | 3 | 0 |

Table 1: Number of regions on the scanned human placenta cryosections positive for MSC-typical marker combinations. All regions listed are CD34- CD45- CD19- CD14- HLA-DR-. One region corresponds to a single cell or to a few cells in very close proximity of each other.

CONCLUSIONS

- To date, it is not known if MSCs express the same set of surface markers in vivo and after in vitro expansion. Here, we present for the first time that CD73+ CD90+ CD105+ CD45- CD34- CD14- CD19- MSCs could be found in human placenta cryosections.
- However, the numbers of triple-positive areas are rather low, suggesting that most MSCs do not express all characteristic markers in vivo but up-regulate their expression during in vitro culture. In line with this, we also detected some additional CD73+ CD90+ and very few CD73+ CD105+ double-positive lineage marker-negative regions, which might be precursors or subsets of MSCs.
- Noteworthy, our method allows different applications to further understand the role and function of MSCs in vivo. For example, it offers the opportunity to track injected MSCs in patient biopsies and thereby understand their fate and interaction with recipient immune cells after application.
- Chipcytometry is useful tool in translational pathology. Further research is needed to validate Chipcytometry for use in clinical trial support.